





ENHANCING ANTERIOR RESTORATIONS - The Role of Modern Matrix Systems

Anterior Matrix System aids in creating precise contours and contacts in restorative procedures for anterior teeth, enhancing both the functional and aesthetic outcomes.

1) Transparent Anterior Matrix System

CLASSIFICATION

2) Non-Transparent Anterior Matrix System



Mylar Strip

Application

Smaller defects, where adjacent tooth has flat contact region

Limitation

Can't contour vast areas, as

it is flexible

Introduced by **Dr. David Clark** in 2007.



Bioclear Matrix

Advantages

- **Excellent Cervical adaptation**
- Preserves Gingival papilla
- Incisal tab for proper orientation

Indications

- All Anterior cases
- · Diastema closure
- · Closure of black triangle

Fusion Anterior Matrix

Unica Anterior Matrix

Advantages

laterals & canines.

Advantages

Incisor

Preserves the ideal anatomical

Developed by **Polydentia** in association with Style Italiano.

Indications

- Class III, IV and V Cavity Restorations
- Direct Composite Veeners
- Shape Modification



curvature in gingivo-incisal and faciallingual orientation.

Ideal for restoration of Maxillary Central

Is used for restoration of smaller anterior

teeth like mandibular central incisor, peg

Introduced by Garrison Dental

Indications

- Class III and IV Cavity Restorations
- Deep Restoration
- Composite Veeners

Transparent Crowns Strip crowns / Cellulose crowns

Advantages

Precontouring easily creates occluso-gingival anatomy and avoids flat embrasures

Application

- Class III and Class IV Composite Restoration
- Diastema closure

Burton-band Anterior Matrix

Advantages

- Complete access to restoration
- Simple installation in single step
- Control over contour and proximal contact
- Consistently seals the gingival margins

Designed by Dr.Mathew Burton.

Appearance

- Metal Matrix is affixed to a plastic wedge.
- Wedge's notched edge clicks into the position and stabilizes the matrix.

Contoured Anterior Matrix

Advantages

Smaller defects, where adjacent tooth has flat contact region

Blue View Vari-Strip

Application

Can't contour vast areas, as it is flexible



References - Urkande NK, Mankar N, Nikhade PP, Chandak M, Ikhar A, Patel A. Anterior Matrix Systems for Composite Restorations: A Review. Cureus. 2023 Apr 4;15(4):e37145. doi: 10.7759/cureus.37145