

Innovations in Dental Matrices: Temporary Contours for Lasting Smiles

The Composi-Tight® 3D Fusion™ sectional matrix system. Designed to tackle complex Class II restorations with precision.



Garrison matrix

The SuperMat is a Universal Matrix Tensioning System designed for large areas



Supermat matrix

All in one Matrix system. Innovative integrated tightening /opening system.



Metafix matrix

- Dental matrix is introduced by **Dr. Louis Jack** in **1871.**

OBJECTIVES OF DENTAL MATRICES

- Assure dryness & non-contamination of the details & the space to be covered with & occupied by the setting restorative material
- Maintain its shape during hardening of the material
- Confine the restorative material within the cavity preparation & predetermined surface configuration



Unica matrix

Ideal matrix for class III,IV restorations . Composite veneers . It allows to reconstruct interproximal and cervical margins with a single matrix



360 degree cervical matrix

The matrix rotates 360 degrees on its handle It is clear for cutting and flexible to facilitate a well sealed gingival margin .



Bioclear matrix

cutting-edge technique for non-invasive anterior, posterior composite restorations, including an innovative approach to black triangle closures.